

All products can be customized., It is now produced for vaccines.
 For orders of 100M per year, the purchase puts the company logo.

Exclusive contracts are also possible. have 510K CRT

The amount is 05.FEB.2021




Quoted price and production capacity:

You can choose between LC at Sight or 10% down payment.

Syringes ML	Needle Gauge	Needle Length	LL		Certification	LS	
			without needle	with Hypodermic needle		without needle	with Hypodermic needle
1ml	25G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9001/9006		●
	23G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9002		●
	25G	1.5"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9003		●
	23G	1.5"		●	(Give us available certification)		●
	///	///		●	GMP-FDA,TCVN,SIO 14001:2015, EN ISO		●
3ml	18G	1.5"		●	FDA/510K		●
	20G	1"		●	FDA/510K		●
	22G	1.5"		●	FDA/510K		●
	25G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9004/9001		●
	23G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9002/9006		●
	25G	1.5"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9003		●
	23G	1.5"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9005		●
	27G	0.5"		●	(Give us available certification)		●
	25G	5/8"		●	FDA/510(K)		●
	///	///		●	GMP-FDA,TCVN,SIO 14001:2015, EN ISO		●
5mL	25G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9006/9004		●
	23G	1"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9007/9005		●
	25G	1.5"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9006		●
	23G	1.5"		●	FDA/CE/ 510(K)/ GMP/WHO/ISO-13485/ ISO-9008		●

Product	Specification	Certification	
Precision Glide needle	Part #302804, .6mm x 25mm	GMP-FDA,TCVN, SIO 14001:2015, EN ISO	
Precision Glide needle	Part #305902, .6mm x 25mm	GMP-FDA,TCVN, SIO 14001:2015, EN ISO	
Oxygen Pipes	////	GMP-FDA,TCVN, SIO 14001:2015, EN ISO	



Specification	Type	Luer Model	Packing
0.5ml	Three Part	Luer Slip	100Pcs*36Boxes
1ml	Three Part	Luer Slip	100Pcs*36Boxes
2ml	Three Part	Luer Slip	100Pcs*24Boxes
2.5ml	Three Part	Luer Slip	100Pcs*24Boxes
3ml	Three Part	Luer Slip	100Pcs*24Boxes
5ml	Three Part	Luer Slip	100Pcs*18Boxes
10ml	Three Part	Luer Slip	100Pcs*12Boxes
20ml	Three Part	Luer Slip	100Pcs*8Boxes
30ml	Three Part	Luer Slip	70Pcs*8Boxes
50ml	Three Part	Luer Slip	45Pcs*8Boxes
60ml	Three Part	Luer Slip	45Pcs*8Boxes
2ml	Three Part	Luer Lock	100Pcs*24Boxes
2.5ml	Three Part	Luer Lock	100Pcs*24Boxes
3ml	Three Part	Luer Lock	100Pcs*24Boxes
5ml	Three Part	Luer Lock	100Pcs*18Boxes
10ml	Three Part	Luer Lock	100Pcs*12Boxes
20ml	Three Part	Luer Lock	100Pcs*8Boxes
50ml	Three Part	Luer Lock	45Pcs*8Boxes
60ml	Three Part	Luer Lock	45Pcs*8Boxes
100ml	Three Part	Luer Lock	25Pcs*8Boxes

Auto-disable Syringe



Description	Specification	Type	Luer Model	Packing
Auto Disable Syringe	0.5ml	Three Part	Luer Slip	100Pcs*36Boxes
Auto Disable Syringe	1 ml	Three Part	Luer Slip	100Pcs*36Boxes
Auto Disable Syringe	2ml	Three Part	Luer Lock	100Pcs*24Boxes
Auto Disable Syringe	2.5ml	Three Part	Luer Lock	100Pcs*24Boxes
Auto Disable Syringe	5ml	Three Part	Luer Lock	100Pcs*18Boxes
Auto Disable Syringe	10ml	Three Part	Luer Lock	100Pcs*12Boxes
Auto Disable Syringe	20ml	Three Part	Luer Lock	100Pcs*8Boxes

Auto-disable safety syringe



1. Check whether the individual pack is intact before use. Forbidden to use if the package is damaged, protective cap or position safety card fallen down.
Peel off individual package to take out the product.



2. Take down the protective cap of needle, drawing a certain amount of liquid, put the needle upwards for air exhaust until no bubbles in the liquid and reach its required dosage.



3. Take down position safety card of plunger.



4. Disinfect the puncture site of human body, then inject.




5. After injection, pull out syringe, pull back the plunger, at same time needle will be retracted to syringe barrel.



6. Plunger will be easily broken after pulling back to the bottom. Throw auto-disable syringe into appointed container.

Description	Specification	Type	Luer Model	Packing
Safety Syringe	1ml	Three Part	Luer Slip	100Pcs*36Boxes
Safety Syringe	2ml	Three Part	Luer Lock	100Pcs*24Boxes
Safety Syringe	2.5ml	Three Part	Luer Lock	100Pcs*24Boxes
Safety Syringe	5ml	Three Part	Luer Lock	100Pcs*18Boxes
Safety Syringe	10ml	Three Part	Luer Lock	100Pcs*12Boxes
Safety Syringe	20ml	Three Part	Luer Lock	100Pcs*8Boxes

► Hypodermic Needles and Safety Needles



Specification		Color Code	Needle Length	Packing
0.3	30G	Light Yellow	3/8"	100Pcs*100Boxes
0.33	29G	Red	1/2"	100Pcs*100Boxes
0.36	30G	Blue Green	1/2"	100Pcs*100Boxes
0.41	27G	Grey	1/2"	100Pcs*100Boxes
0.45	26G	Brown	5/8"	100Pcs*100Boxes
0.51	25G	Orange	3/4"	100Pcs*100Boxes
0.56	24G	Purple	3/4"	100Pcs*100Boxes
0.63	23G	Blue	1"	100Pcs*100Boxes
0.71	22G	Blanck	1 1/4"	100Pcs*100Boxes
0.82	21G	Green	1 1/2"	100Pcs*100Boxes
0.9	20G	Yellow	1 1/2"	100Pcs*100Boxes
1.08	19G	Milk White	1 1/2"	100Pcs*100Boxes
1.28	18G	Pink	1 1/2"	100Pcs*100Boxes
1.63	16G	White	1 1/2"	100Pcs*100Boxes