

High use po

CARDINAL HEALTH

Esteem™ Comfort

Esteem™ protection with improved elongation

Cardinal Health Esteem™ Comfort is our newest generation of Esteem™ stretchy nitrile exam gloves. This glove is the perfect combination of improved elongation with outstanding protection from chemotherapy drugs and an FDA approved low dermatitis potential claim.

Features and benefits:

- 30% higher min. elongation¹ than the minimum ASTM standard² for a nitrile exam glove which can help to **reduce cuff tears during donning.**
- Reduces the risk of developing a Type IV Hypersensitivity to Chemical Accelerants.
- Tested against chemotherapy drugs³ (with a minimum breakthrough detection time for Thiopeta >30mins) (Note: test results are printed directly on the packaging.)
- Tested against 18 chemicals⁴
- Textured fingertips for enhanced grip
- Premium box with dispensing film

High Use: A portfolio of non-sterile exam gloves designed to serve a broad range of everyday



Cardinal Health™ Esteem™ Comfort



- Powder-free
- Not made with natural rubber latex
- Ambidextrous
- Color: Blue
- Nominal fingertip thickness: 4.3mil (0.11mm)
- Textured fingertips for enhanced grip

Cat. No.	Size	Pkg.
88ES01XS	X-Small	200/bx, 10bx/cs
88ES02S	Small	200/bx, 10bx/cs
88ES03M	Medium	200/bx, 10bx/cs
88ES04L	Large	200/bx, 10bx/cs
88ES05XL	X-Large	180/bx, 10bx/cs

For more information or to order, contact your Cardinal Health sales representative. Call 877.CARDINAL or visit cardinalhealth.com.

¹ Elongation (EN): The amount of extension (increase) of an object under stress, usually expressed as a percentage of its original length.

² Ultimate Elongation Before Aging (UEBA) per ASTM D6319

³ This glove has been tested for resistance to permeation of various chemotherapy drugs per ASTM D6319, "Standard Practice for Assessment of Resistance of Gloves to Permeation by